



ELECTRIC HEATING INCUBATOR



BXP series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

ELECTRIC HEATING INCUBATOR >>

BXP series

Product Features >

- LCD display, real-time display of working status;
- BRIGHT I control system, automatically compensate parameters based on the environment;
- With countdown reservation function, the reservation time range is 99 hours 59 minutes;
- Programmable design: 10 segments 100 cycles program setting;
- Double over temperature alarm function (A. controlled by temperature sensor system, B. independent temperature control switch);
- The 6-stage speed regulation of the fan,(BXP-N series no fan);
- Preheating chamber design, to ensure rapid heating and good heat distribution effect;
- Double door design, with glass inner door;
- Equipped with 1 RS-485 interface.



Lock



Fan



Shelf

Model	BXP-65	BXP-130	BXP-210	BXP-280	BXP-450
Volume [L]	64	130	210	280	452
Temp. Range [°C]	RT+5~100				
Temp. Fluctuation at 37°C [°C]	±0.3				
Temp. Uniformity at 37°C [°C]	±0.8				
Heating-up Time to 60°C [min]	16	20	15	15	25
Recovery time after 30 seconds door open at 60°C [min]	14	13	10	10	11
Power [W]	500	500	750	950	1500
Internal Dimensions W×D×H [mm]	410×390×450	570×400×570	630×530×630	720×530×745	600×580×1300
Housing Dimensions W×D×H [mm]	575×566×722	715×616×842	775×746×902	865×730×1017	774×845×1704
Number of Shelves (std./max.) [pc]	2/5	2/7	3/7	3/9	4/13
Net Weight /Gross Weight [kg]	46/64	57/75	77/96	97/118	151/180
Package Size [mm]	675*726*872	835*776*992	895*906*1052	985*906*1167	894*1005*1854