

HORIZONTAL SHAKER



BXYC-WX-L Series



BXYC-WX-C Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

HORIZONTAL SHAKER >>

BXYC-WX-L Series

Product Feature

- The characteristics of the illuminated horizontal intelligent precision shaking table product are the same as those of BXYC-WX series;
- Adopting full spectrum high intensity LED lights, with infinitely adjustable lighting intensity.



LED Light

BXYC-WX-C Series

Product Feature

- The characteristics of the carbon dioxide horizontal intelligent precision shaking table are the same as those of the BXYC-WX series;
- Equipped with high-precision carbon dioxide infrared sensor and pressure reducing valve.



Model	BXYC-WX2350L/2280L/2200L BXYC-WX2350C/2280C/2200C	BXYC-WX1350L/1280L BXYC-WX1350C/1280C
Temp. Range [°C]	4~60	4~60
Temp. Resolution [°C]	0.1	0.1
Temp. Fluctuation at 37°C [°C]	±0.3	±0.3
Temp. Uniformity at 37°C [°C]	±1	±1
Humidity Range [%]	-	Ambient Humidity-95%RH
Humidity Fluctuation [%]	-	±3
CO2 Range [%]	—/0~20	—/0~20
CO2 Accuracy [%]	—/0.1	—/0.1
Illuminance [Lux]	0~20000/—	0~20000/—
Rotation Rate [rpm]	30~300	30~300
Rotation Accuracy [rpm]	±1	±1
Vibration Amplitude [mm]	∅26	∅26
Timer Range [h]	0~999	0~999
Power Supply [V, Hz]	220, 50/60	220, 50/60
Power [W]	1800/1600/1500	2800/2400
Internal Dimensions W×D×H [mm]	1035×675×470/1010×590×470/1010×590×310	1035×675×470/1010×590×470
Housing Dimensions W×D×H [mm]	1380×800×1010/1350×700×1000/1260×730×850	1380×800×1010/1350×700×1000
Tray Size W×D [mm]	970×605/950×530/950×530	970×605/950×530
Standard Configuration [ml×pc]	2000×8/2000×6/500×20	2000×8/2000×6
Full Load Configuration [ml×pc]	600	50×103/100×103/250×63/500×41/1000×22/2000×15/3000×11/5000×7
	400	50×71/100×71/250×45/500×37/1000×18/2000×15 or 3000×8 or 5000×7
	200	50×71 or 100×71 or 250×45 or 500×37 or 1000×18
Net weight [kg]	278/262/228	282/266
Gross weight [kg]	358/322/256	362/328
Package size [mm]	1540*1000*1240/1510*910*1190/1400*900*1030	1540*1000*1240/1510*910*1190